

Antibody Identification

Session 2/3 – Putting it into practice

Case panels

Session 2/3 – Putting it into practice
Case 1 Panel

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	ENZ
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	3	0	+	0	0	0	+	0	+	0	+		0	0
2	R ₁ R ₁	+	+	0	0	+	+	+	0	+	1	0	0	+	0	+	0	+	0	+	0	HLA+	1	0
3	R ₂ R ₂	+	0	+	+	0	0	+	0	+	3	0	0	+	0	+	0	0	+	+	0		2	0
4	r'r	0	+	0	+	+	0	+	0	+	0	0	0	+	0	+	0	+	0	0	+		2	0
5	r''r	0	0	+	+	+	+	0	+	0	1	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	0	+	0	+	4	0	+	0	0	0	+	0	+	+	0	HLA+	2	0
7	rr	0	0	0	+	+	+	+	0	+	0	0	+	+	0	0	+	W	0	+	0		1	0
8	rr	0	0	0	+	+	+	0	0	+	0	0	0	+	+	0	+	+	0	0	+		0	0
9	rr	0	0	0	+	+	0	+	+	0	1	0	0	+	0	0	+	+	0	+	0		1	0
10	rr	0	0	0	+	+	+	0	0	+	3	+	0	+	0	0	+	0	+	+	0		0	0
																						Auto	0	/
																						K control	2	/

Session 2/3 – Putting it into practice
Case 2 Panel

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0		0	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		4	0
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		4	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		4	0
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		0	0
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		0	0
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		4	0
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		4	0
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		0	0
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		0	0
																						Auto	0	/
																						K control	2	/

Session 2/3 – Putting it into practice
Case 3 Panel

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	EIAT
1	$R_1^w R_1$	+	+	0	0	+	+	0	+	+	2	0	0	+	0	0	+	0	+	+	0		0	0
2	$R_1 R_1$	+	+	0	0	+	0	+	0	+	0	+	+	+	0	+	0	+	0	0	+		0	0
3	$R_2 R_2$	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	0		4	5
4	$r'r$	0	+	0	+	+	0	+	0	+	0	0	0	+	0	0	+	0	+	+	0		2	5
5	$r''r$	0	0	+	+	+	0	+	+	0	4	0	0	+	0	0	+	+	0	0	+		4	5
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		3	5
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	+	+	+	0		3	5
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		3	5
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		3	5
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		3	5
																						Auto	0	/
																						K control	2	/

Session 2/3 – Putting it into practice
Case 4 Panel

Cell	Rh	D	C	E	c	e	M	N	S	s	P1	Lu ^a	K	k	Kp ^a	Le ^a	Le ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Other	IAT	EIAT
1	R ₁ ^w R ₁	+	+	0	0	+	+	0	+	0	0	0	0	+	0	0	+	+	0	+	0		0	0
2	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	0	+	0	+		0	0
3	R ₂ R ₂	+	0	+	+	0	+	0	+	0	0	0	0	+	0	+	0	0	+	+	+		0	0
4	r'r	0	+	0	+	+	0	+	+	0	0	0	0	+	0	0	+	0	+	+	0		2	4
5	r''r	0	0	+	+	+	0	+	0	+	4	0	0	+	0	0	+	+	0	0	+		2	4
6	rr	0	0	0	+	+	+	0	+	0	4	0	+	0	0	0	+	+	0	0	+		2	4
7	rr	0	0	0	+	+	0	+	0	+	2	+	+	+	0	+	0	0	+	+	0		2	4
8	rr	0	0	0	+	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+		2	4
9	rr	0	0	0	+	+	+	0	0	+	2	0	0	+	0	+	0	+	0	0	+		2	4
10	rr	0	0	0	+	+	0	+	0	+	3	0	0	+	0	+	0	+	0	+	0		2	4
																						Auto	0	/
																						+ control	2	/